



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Weighting wire, 30 cm long, 114 g

€55.00
plus VAT &
Shipping

Product Images



Description

Flexible lead wire with resistant vinyl coating,
can be formed into any desired shape.
Particularly suitable for weighing down sample vessels in a water bath.

Wire thickness: 6.4 mm

Length: 30 cm

Weight: approx. 114 g

Additional Information

No.	6-2335
Manufacturer (Brand)	neoLab
EAN	4058072025748
Transport temperature	Room temperature
Color	red
Material	Lead
Form	annular. open
Length	300 mm
Weighting material properties	with special chemical resistance
Weight	114 g
Thickness	6.4 mm

