



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Lead ring with vinyl coating, 51 mm inner Ø, for 250 to 1000 ml

€40.10
plus VAT &
Shipping

Product Images



Description

Lead rings for weighting glass vessels in a water bath.
To reduce glass breakage or leakage of chemicals.
With vinyl coating, particularly resistant to chemicals.

Additional Information

No.	6-1012
Manufacturer (Brand)	neoLab
EAN	4058072024642
Transport temperature	Room temperature
Color	blue
Material	Lead Polyvinyl chloride (PVC. unspecified)
Form	annular. closed
Weighting material properties	with special chemical resistance
DM inside	51 mm
Weight	454 g

