



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® lead ring with vinyl coating, 70 mm inner Ø, for 1 to 4 l

**€51.90**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Lead rings for weighting glass vessels in a water bath.  
To reduce glass breakage or leakage of chemicals.  
With vinyl coating, particularly resistant to chemicals.

## Additional Information

---

No.	6-1014
Manufacturer (Brand)	neoLab
EAN	4058072024666
Transport temperature	Room temperature
Color	blue
Material	Lead Polyvinyl chloride (PVC. unspecified)
Form	annular. closed
Weighting material properties	with special chemical resistance
DM inside	70 mm
Weight	908 g

