



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Weighting ring 1200 g, 100 mm inner diameter.

€90.00
plus VAT &
Shipping

Product Images



Description

For weighing down glass vessels to prevent them from falling over or floating up in the water bath. The rings made of lead with PVC coating are simply slipped over the neck for weighting.

Additional Information

No.	3-1630
Manufacturer (Brand)	neoLab
EAN	4058072013240
Transport temperature	Room temperature
Color	blue
Material	Lead Polyvinyl chloride (PVC. unspecified)
Form	annular. open
DM inside	100 mm
Weight	1200 g

