



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



SERVA PBS buffer (10x) sterile 1 l

€76.00
plus VAT &
Shipping

Product Images



Description

10x concentrated phosphate buffered salt solution, sterile filtered. PBS buffer is a commonly used buffer in protein detection systems, such as Western blot, ELISA and other enzyme assays, for immunocytological and immunohistochemical detection, in situ hybridizations, apoptosis assays and nuclear staining. 1x PBS is often used as a protein solvent and diluent.

pH: 7.2-7.6

Additional Information

| | |
|-----------------------|------------------|
| No. | SV-0092 |
| Manufacturer (Brand) | SERVA |
| EAN | 4058072372293 |
| Transport temperature | Room temperature |

