



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## SERVA VisiBlot Standard I, 500 µl

**€283.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

VisiBlot Standard I is a mixture of recombinant proteins in the molecular weight range from 25 kDa to 150 kDa. The protein bands at 25 kDa, 45 kDa and 85 kDa are pre-stained for better visibility during electrophoresis. The other 5 proteins possess IgG binding sites, so these proteins are easily detectable via antibody detection on the Western blot membrane. Since these proteins are not covalently labeled with dyes, they allow accurate assignment of molecular weight. The recommended application volume for a minigel is 5 µl/trace.

- Ready-to-use, no reconstitution, dilution or heating required
- Pre-stained bands for better visibility during electrophoresis and on the blot membrane
- Detection of marker proteins on the blot membrane according to the detection system of the antibodies, e.g. horseradish peroxidase, alkaline phosphatase
- Simple molecular weight determination of the detected protein on the blot membrane

Molecular weight range (kDa): 25-150 with IgG binding sites (plus 3 colorimetric bands)

## Additional Information

---

No.	SV-0055
Manufacturer (Brand)	SERVA
EAN	4058072371920
Transport temperature	Room temperature

