



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine insert for 4 x 4 reaction tubes 0.5 ml, acrylic transparent

€34.10
plus VAT &
Shipping

Product Images



Description

Insert racks made of acrylic glass for reaction tubes with 0.5 ml content.
Dimensions: approx. 75 x 75 x 45 mm (L x W x H)

Additional Information

No.	2-5019
Manufacturer (Brand)	neoLabLine
EAN	4058072003845
Transport temperature	Room temperature
Color	transparent (clear)
Material	Polymethylmetacrylate (PMMA) (acrylic glass)
Wide	75 mm
Height	45 mm
Length	75 mm
Depth	75 mm
suitable for	Micro Vessels / Tubes
Squares (WxL)	4 x 4
for vessels with DM up to	5.5 mm
for vessels with volume up to	0.5 ml
for type (micro) vessel	Reaction vessels

