



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® desiccator insert made of stainless steel, material: 1.4301, type 304, rustproof

**€27.25
plus VAT &
Shipping**

Product Images



Description

Stainless steel desiccator insert, DN 200 The desiccator insert is a sieve plate that separates the desiccator base into two sections - a drying agent is placed in the lower part, and the substance to be dried is placed in the upper part. The substance to be dried is placed on the sieve plate in an appropriate vessel for drying. Special features and advantages: The heat-resistant material allows hot samples to be placed in the desiccator Made of heat-resistant material Made of stainless steel Application: Drying of materials

Additional Information

No.	SD-2607
Manufacturer (Brand)	DURAN
VGKL number	122908061
EAN	4058072362416
Manufacturer EAN	4032051061651
Gross weight	0.146
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.14600 kg
diameter	190 mm
Material	Stainless steel
Country of origin	Germany
sterile	No
Region of origin	Hesse
Customs tariff number	72193410
ND (nominal diameter)	DN 200

