



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® desiccator insert made of stainless steel, material: 1.4301, type 304, rustproof

**€20.18
plus VAT &
Shipping**

Product Images

Description

Stainless steel desiccator insert, DN 150 The desiccator insert is a sieve plate that separates the desiccator base into two sections – a drying agent is placed in the lower part, and the substance to be dried is placed in the upper part. The substance to be dried is placed on the sieve plate in an appropriate vessel for drying. Special features and advantages: The heat-resistant material allows hot samples to be placed in the desiccator Made of heat-resistant material Made of stainless steel Application: Drying of materials

Additional Information

No.	SD-2606
Manufacturer (Brand)	DURAN
VGKL number	122908057
EAN	4058072362409
Manufacturer EAN	4032051061644
Gross weight	0.08
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.08000 kg
diameter	140 mm
Material	Stainless steel
Country of origin	Germany
sterile	No
Region of origin	Hesse
Customs tariff number	76169950
ND (nominal diameter)	DN 150

