



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® desiccator insert made of porcelain

€211.53
plus VAT &
Shipping

Product Images



Description

Desiccator insert made of porcelain, DN 300 The desiccator insert is a sieve plate that separates the desiccator base into two parts - a drying agent is placed in the lower part, and the substance to be dried is placed in the upper part. The substance to be dried is placed on the sieve plate in an appropriate vessel for drying. Special features and advantages: The heat-resistant material allows hot samples to be placed in the desiccator Made of heat-resistant material Application: drying of materials

Additional Information

No.	SD-2835
Reg. trademark	DURAN ®
Manufacturer (Brand)	DURAN
VGKL number	122972569
EAN	4058072364281
Manufacturer EAN	4032051064768
Gross weight	1.358
Transport temperature	Room temperature
Storage temperature	Room temperature
Net weight	1.33800 kg
diameter	280 mm
Material	Porcelain
sterile	No
Country of origin	Germany
Packaging volume	0.00078400 ccm
Customs tariff number	69099000
ND (nominal diameter)	DN 300

