



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® desiccator insert made of porcelain

€112.41
plus VAT &
Shipping

Product Images

Description

Desiccator insert made of porcelain, DN 200 The desiccator insert is a sieve plate that separates the desiccator base into two parts - a drying agent is placed in the lower part, and the substance to be dried is placed in the upper part. The substance to be dried is placed on the sieve plate in an appropriate vessel for drying. Special features and advantages: The heat-resistant material allows hot samples to be placed in the desiccator Made of heat-resistant material Application: drying of materials

Additional Information

No.	SD-2833
Manufacturer (Brand)	DURAN
VGKL number	122972561
EAN	4058072364267
Manufacturer EAN	4032051064744
Gross weight	0.528
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.50900 kg
diameter	184 mm
Material	Porcelain
Country of origin	Czech Republic
sterile	No
Packaging volume	0.00036100 cc
Customs tariff number	69091100
ND (nominal diameter)	DN 200

