



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## DURAN® desiccator insert made of porcelain

**€37.25**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Desiccator insert made of porcelain, DN 100 The desiccator insert is a sieve plate that separates the desiccator base into two parts - a drying agent is placed in the lower part, and the substance to be dried is placed in the upper part. The substance to be dried is placed on the sieve plate in an appropriate vessel for drying. Special features and advantages: The heat-resistant material allows hot samples to be placed in the desiccator Made of heat-resistant material Application: drying of materials

## Additional Information

---

No.	SD-2831
Manufacturer (Brand)	DURAN
VGKL number	122972546
EAN	4058072364243
Manufacturer EAN	4032051064720
Gross weight	0.06
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.05700 kg
diameter	90 mm
Material	Porcelain
Country of origin	Czech Republic
sterile	No
Customs tariff number	69091100
ND (nominal diameter)	DN 100

