



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® micro inserts clear glass 0.1 ml,
conical, 13 mm tip, 100 pcs/pack**

**€26.40
plus VAT &
Shipping**



Description

Conical microinsert

13 mm tip

Microinserts can be fixed in the transparent PE stoppers: in this way, no springs are required.

Suitable for flat-bottomed vials Order No. 7-0739 and 7-0740

Additional Information

| | |
|-----------------------|--|
| No. | 7-0741 |
| Manufacturer (Brand) | neochrom |
| EAN | 4058072029944 |
| Transport temperature | Room temperature |
| diameter | 5 mm |
| Material | Borosilicate glass (1st hydrolytic class) |
| Volume | 0.15 ml |
| Height | 34 mm |
| Floor outside | conical. pointed |
| Length tip | 13 mm |
| for type vial/vial | for flat-bottomed vials 1 ml with insertion lock |

