



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neochrom® micro insert clear glass, for thread fl. ND13, conical, 15 mm tip, 100 pcs/pack

Special Price

€14.00 was

€21.00

**plus VAT &
Shipping**



Description

Conical micro insert for threaded vials ND 13
15 mm tip,
capacity 0.3 ml, 40 x 6 mm Ø.

Additional Information

No.	7-0768
Manufacturer (Brand)	neochrom
EAN	4058072030179
Transport temperature	Room temperature
diameter	6 mm
Material	Borosilicate glass (1st hydrolytic class)
Volume	0.4 ml
Height	40 mm
Floor outside	conical. pointed
Length tip	15 mm
for type vial/vial	for threaded vials. 13 mm

