



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Cryomaster® perforated stainless steel plate 10x10 for stainless steel cryobox 136x136x53 mm

**Price on
request
plus VAT &
Shipping**

Product Images

Description

Cryomaster® perforated plate made of stainless steel 10x10

The combination of Cryomaster cryobox and perforated plate made of high-quality stainless steel is characterized by its robustness in contrast to cryoboxes made of cardboard or plastic. This ensures a long-lasting product and is therefore ideally suited for long-term storage of important samples.

- for 136 x 136 mm cryoboxes

Additional Information

No.	4-6025
Manufacturer (Brand)	Cryomaster
EAN	4058072191382
Transport temperature	Room temperature
Color	silver
Material	Stainless steel
Wide	134 mm
Length	134 mm
Depth	134 mm
Squares (WxL)	10 x 10

