



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e



Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine perforated plate transparent, coded, grid 9x9 for neoBox 81 plus

€1.51
plus VAT &
Shipping

Product Images



Description

Numbered, coded perforated plate, suitable for neoBox 81 plus.

The perforated plate is inserted into the lower part, so that each "specimen position" is clearly numbered and assigned. Thus, all 81 positions (grid 9 x 9) in the box are unambiguous and there is no mix-up of samples.

The numbered perforated plate is made of transparent polycarbonate and is temperature resistant from -90°C to +121°C and autoclavable.

Additional Information

No.	2-2939
Manufacturer (Brand)	neoLabLine
EAN	4058072168100
Transport temperature	Room temperature
autoclavable at 121°C	No
Material	Polycarbonate (PC)
Temp MAX	121 °C
Temp MIN	-90 °C
sterile	No
Inserts properties	alphanumeric coded
Squares (WxL)	9 x 9

