



WTW CelloX® 325

Product Images



Description

Galvanic oxygen sensor with OxiCal®-SL calibration and storage vessel, waterproof plug (IP 67), cable length 1.5 m. Incl. cleaning solution, grinding foil, 3 interchangeable membrane heads and electrolyte.

Additional Information

No.	HT-0209
Reg. trademark	OxiCal ®
Manufacturer (Brand)	WTW
VGKL number	688132533
EAN	4066292544321
Gross weight	0.55
Preference identifier	Ja
UN number	2790
Country of origin	Deutschland
Customs tariff number	90279000

