



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



SI Analytics ELECTRODE PH COMBINATION N 6280 120MM

€259.00
plus VAT &
Shipping

Product Images

Description

pH combination electrode with plug-in head N 6280, glass shaft, ceramic diaphragm, electrolyte KCl 3 mol/l, Silamid® derivative, spherical membrane, A-glass, length 120 mm, 12 mm Ø, -5...+100 °C, 0...14 pH

Additional Information

No.	SY-0158
Reg. trademark	Silamide ®
Manufacturer (Brand)	SI Analytics
VGKL number	612720015
EAN	4058072172732
Gross weight	0.105
Transport temperature	Room temperature
Net weight	0.10500 kg
Country of origin	Germany
Customs tariff number	90279000
Length	120 mm
DM Electrode/ Sensor	12 mm
pH MAX	14.0 pH
pH MIN	0.0 pH
Material diaphragm	Ceramics
Material electrode / sensor	Glass
Type plug	Plug-in head

