



WTW F 800 DIN

Product Images

Description

Fluoride selective combination electrode with 1 m DIN cable, 0.02 mg/l to saturation, requires ELY/BR/503 reference electrolyte, ionic strength adapting solution TISAB, standard solution ES/F. Electrolyte included in delivery!

Additional Information

No.	HT-0074
Manufacturer (Brand)	WTW
VGKL number	108116667
EAN	4058072981532
Gross weight	0.28
Country of origin	USA
Packing width	0.23000 m
Packing height	0.06000 m
Packing depth	0.09500 m
Packaging volume	0.00131100 ccm
Customs tariff number	90279000

