



Mettler-Toledo InLab Easy pH electrode

Product Images



Description

Low cost combined pH electrode. Ideal for field measurements and educational use
Specifications: Rugged polysulfone shaft (12 mm diam./ 120 mm long) S7 screw head Ceramic diaphragm Gel electrolyte Ag/AgCl- Reference pH range: 0...14 Temp. range: 0...80°C

Additional Information

No.	TD-0410
Manufacturer (Brand)	Mettler-Toledo
VGKL number	105003010
EAN	4066292343504
Gross weight	0.068
Country of origin	Schweiz
Customs tariff number	90278020

