



Mettler-Toledo InLab Micro pH electrode

Product Images



Description

Designed for measurements in very small vessels such as vials or cuvettes. Specifications: Glass shaft (3 mm diam./ 60 mm long) S7 screw head Ceramic diaphragm ARGENTHAL™- reference with silver ion barrier pH range: 0...14 Temp. range: 0...100°C Incl. 25 ml reference electrolyte (3 mol/L KCl)

Additional Information

No.	TD-0431
Manufacturer (Brand)	Mettler-Toledo
VGKL number	105003160
EAN	4066292343689
Gross weight	0.11
Country of origin	Schweiz
Customs tariff number	90278020

