



Mettler-Toledo InLab Reference

Product Images



Description

Reference electrode for measurements with pH and ion selective half cells.Specifications: Glass shaft (12 mm diam./ 120 mm long)
S7 screw head ARGENTHAL™- reference with silver ion barrier Kerami diaphragm Temp. range: 0...100°C

Additional Information

No.	TD-0442
Manufacturer (Brand)	Mettler-Toledo
EAN	4066292343795
Gross weight	0.15
Country of origin	Schweiz
Customs tariff number	90278020

