



Mettler-Toledo InLab Science pH electrode

Product Images



Description

Ideal for slurries and emulsions
Specifications: Glass shaft (12 mm diam./ 120 mm long) S7 screw head Movable ground-joint diaphragm ARGENTHAL™- reference with silver ion barrier pH range: 0...12 Temp. range: 0...100°C Incl. 25 ml reference electrolyte (3 mol/L KCl)

Additional Information

No.	TD-0418
Manufacturer (Brand)	Mettler-Toledo
VGKL number	105003070
EAN	4058072982591
Gross weight	0.152
Country of origin	Schweiz
Customs tariff number	90278020

