



Mettler-Toledo InLab® Expert Pro-ISM

Product Images



Description

pH electrode InLab® Expert Pro-ISM for SevenCompact instruments: - robust 3in1 electrode with PEEK shaft - polymer electrolyte - maintenance-free - contamination-resistant due to open connection - integrated temperature sensor - Intelligent Sensor Management (ISM): sensor data, calibration, and maximum temperature are stored in the sensor and are automatically transferred to the meter when connected - fixed cable: 1.5m - connector: BNC, Cinch/RCA

Additional Information

No.	TD-0018
Reg. trademark	InLab ®
Manufacturer (Brand)	Mettler-Toledo
VGKL number	105004096
EAN	4066292339941
Gross weight	0.18
Country of origin	Schweiz
Customs tariff number	90278020

