



Mettler-Toledo InLab® Micro Pro-ISM

Product Images



Description

InLab® Micro Pro-ISM - combined pH-micro-electrode with glass shaft and integrated temperature sensor - ideal for small sample volumes / vessels - shaft diameter 5mm - MultiPin™ plug-in head - ARGENTHAL™ - reference system (silver ion free) - Intelligent Sensor Management (ISM): Sensor data, calibration, and maximum temperature are stored in the sensor and are automatically transferred to the meter when connected - Ideal for SevenGo and SevenCompact instruments, necessary connection cable: # 51344291 Technical data - pH range 0 ... 14 - Temperature range °C 0 ... 100 - Temperature sensor NTC 30 kOhm - Reference electrolyte 3 mol/L KCl (#51350072) - Shaft length 130 mm - Shaft diameter 5 mm

Additional Information

No.	TD-0479
Reg. trademark	ARGENTHAL ®
Manufacturer (Brand)	Mettler-Toledo
VGKL number	105004163
EAN	4058072983154
Gross weight	0.14
Country of origin	Schweiz
Customs tariff number	90278020

