



Mettler-Toledo InLab® Pure Pro-ISM

Product Images



Description

InLab® Pure Pro-ISM - combined pH electrode with glass shaft ideal for low conductive media - fixed ground-joint diaphragm - exchangeable bridge electrolyte - integrated temperature sensor - MultiPin™ plug-in head - ARGENTHAL™ - reference system with silver ion barrier - Intelligent Sensor Management (ISM): Sensor data, calibration, and maximum temperature are stored in the sensor and are automatically transferred to the meter when connected - Ideal for SevenGo and SevenCompact instruments, necessary connection cable: # 51344291 Technical data - pH range 1 ... 11 - Temperature range °C 0 ... 80 - Temperature sensor NTC 30 kOhm - Membrane resistance (25 °C) < 50 MOhm - Reference electrolyte 3 mol/L KCl (#51350072) - Bridge electrolyte 1 mol/L KCl (#51343181) - Shaft length 120 mm - Shaft diameter 12 mm

Additional Information

No.	TD-0480
Reg. trademark	ARGENTHAL ®
Manufacturer (Brand)	Mettler-Toledo
VGKL number	105004172
EAN	4066292344143
Gross weight	0.17
Country of origin	Schweiz
Customs tariff number	90278020

