



Mettler-Toledo InLab® Science Pro-ISM

Product Images



Description

InLab® Science Pro-ISM - combined pH electrode with glass shaft and ground-joint diaphragm - integrated temperature sensor - MultiPin™ plug-in head - ARGENTHAL™ reference system with silver ion barrier - Intelligent Sensor Management (ISM): sensor data, calibration, and maximum temperature are stored in the sensor and are automatically transferred to the meter when connected - ideal for SevenGo and SevenCompact instruments, necessary connection cable: # 51344291 Technical data - pH range 0 ... 12 - Temperature range °C 0 ... 100 - Temperature sensor NTC 30 kOhm - Diaphragm resistance (25 °C) < 600 MOhm - Reference electrolyte 3 mol/L KCl (#51350072) - Shaft length 170 mm - Shaft diameter 12 mm

Additional Information

No.	TD-0464
Reg. trademark	ARGENTHAL ®
Manufacturer (Brand)	Mettler-Toledo
VGKL number	105004072
EAN	4066292344006
Gross weight	0.16
Country of origin	Schweiz
Customs tariff number	90278020

