



Mettler-Toledo LE 410, pH electrode

Product Images



Description

3-in-1 pH glass electrode with integrated temperature probe and liquid electrolyte,for easy and fast measurements in the laboratory, aqueous samples like acids, bases etc., electrolyte refillable, fixed cable with BNC/Cinch connector (Available from Q2'2012)

Additional Information

No.	TD-0626
Manufacturer (Brand)	Mettler-Toledo
VGKL number	105000348
EAN	4066292345492
Gross weight	0.17
Country of origin	China
Customs tariff number	90278020

