



WTW SenTix® 22

Product Images

Description

Plastic pH electrode with gel electrolyte and BNC connector, without temperature probe, 1 m cable

Additional Information

No.	HT-0007
Manufacturer (Brand)	WTW
VGKL number	618113633
EAN	4066292542570
Gross weight	0.14
Preference identifier	Ja
Country of origin	Deutschland
Packing width	0.25500 m
Packing height	0.05500 m
Packing depth	0.05500 m
Packaging volume	0.00077138 ccm
Customs tariff number	90279000

