



WTW SenTix® 41

Product Images



Description

Plastic pH electrode with gel electrolyte, integrated temperature sensor, waterproof DIN plug and banana plug, 1 m cable

Additional Information

No.	HT-0008
Manufacturer (Brand)	WTW
VGKL number	618121635
EAN	4058072171230
Gross weight	0.18
Preference identifier	Ja
Country of origin	Deutschland
Packing width	0.25500 m
Packing height	0.05500 m
Packing depth	0.05500 m
Packaging volume	0.00077138 ccm
Customs tariff number	90279000

