



WTW SenTix® 51

Product Images

Description

Epoxy pH electrode with liquid electrolyte, integrated temperature sensor, waterproof DIN plug and banana plug, 1 m cable

Additional Information

No.	HT-0022
Manufacturer (Brand)	WTW
VGKL number	618113651
EAN	4066292542693
Gross weight	0.21
Preference identifier	Ja
Country of origin	Deutschland
Packing width	0.25500 m
Packing height	0.05500 m
Packing depth	0.05500 m
Customs tariff number	90279000

