



WTW SenTix® 60

Product Images



Description

Glass pH electrode with liquid electrolyte, S7 plug-in head, length 120 mm. Without connection cable.

Additional Information

No.	HT-0011
Manufacturer (Brand)	WTW
VGKL number	618113639
EAN	4058072177218
Gross weight	0.15
Preference identifier	Ja
Country of origin	Deutschland
Packing width	0.25500 m
Packing height	0.05500 m
Packing depth	0.05500 m
Packaging volume	0.00077138 ccm
Customs tariff number	90279000

