



WTW SenTix® 81

Product Images



Description

Glass pH electrode, liquid electrolyte, integrated temperature sensor, waterproof DIN plug and banana plug, length 120 mm, 1 m cable

Additional Information

No.	HT-0014
Manufacturer (Brand)	WTW
VGKL number	618141014
EAN	4058072187736
Gross weight	0.21
Preference identifier	Ja
Country of origin	Deutschland
Packing width	0.25500 m
Packing height	0.05500 m
Packing depth	0.05500 m
Packaging volume	0.00077138 ccm
Customs tariff number	90279000

