



WTW SenTix® 980

Product Images



Description

IDS pH electrode with glass shaft, liquid electrolyte (3 mol/l KCl), platinum diaphragm, built-in temperature sensor, QSC-capable, memory for 10 calibrations, stored serial number, 1.5 m cable with digital plug for MultiLine® and inoLab® IDS.

Additional Information

No.	HT-0052
Manufacturer (Brand)	WTW
VGKL number	108129800
EAN	4066292542976
Gross weight	0.2
Country of origin	Deutschland
Packing width	0.28000 m
Packing height	0.06500 m
Packing depth	0.10500 m
Packaging volume	0.00191100 ccm
Customs tariff number	90278020

