



WTW SenTix® 980-P

Product Images

Description

IDS pH electrode with glass shaft, liquid electrolyte (3 mol/l KCl), platinum diaphragm, built-in temperature sensor, QSC-capable, memory for 10 calibrations, stored serial number, plug-in head for AS IDS/1.5 cable or radio module.

Additional Information

No.	HT-0046
Manufacturer (Brand)	WTW
VGKL number	108123762
EAN	4066292542914
Gross weight	0.27
Country of origin	Deutschland
Customs tariff number	90278020

