



WTW SenTix® HW

Product Images



Description

Precision pH electrode with ground-joint diaphragm for low-ion media and S7 plug-in head, length 165 mm. Without connection cable.

Additional Information

No.	HT-0021
Manufacturer (Brand)	WTW
VGKL number	618113650
EAN	4066292542686
Gross weight	0.16
Preference identifier	Ja
Country of origin	Deutschland
Packing width	0.25500 m
Packing height	0.05500 m
Packing depth	0.05500 m
Packaging volume	0.00077138 ccm
Customs tariff number	90279000

