



WTW SenTix® HW-T 900

Product Images

Description

IDS precision pH electrode with ground-joint diaphragm for low-ion media, length 165 mm, with built-in temperature sensor, QSC-capable, memory for 10 calibrations, stored serial number, 1.5 m cable with waterproof digital plug for MultiLine®, Mono

Additional Information

No.	HT-0043
Manufacturer (Brand)	WTW
VGKL number	108113753
EAN	4066292542884
Gross weight	0.37
Country of origin	Deutschland
Customs tariff number	90278020

