



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® pH electrode Standard T BNC

Special Price

€143.00 was

€183.00

**plus VAT &
Shipping**

Product Images



Description

The renowned Swiss manufacturer has developed a complete line of pH electrodes for glass-bodied laboratories. The new XS sensors are available with S7 head or cable fixation with BNC or DIN connector.

A special polymer electrolyte/KCl or KCl electrolyte solution combined with a high performance reference system produces a full range of electrodes.

Blue inner buffer provides a clear view of the contact with the inner pH glass membrane.

All XS sensors can be connected to commercially available pH meters.

With temperature sensor NTC 30 $K\Omega$

Additional Information

No.	4-0938
Manufacturer (Brand)	neoLab
EAN	4058072189549
Transport temperature	Room temperature
Temp MAX	100 °C
Temp MIN	-10 °C
Length	120 mm
suitable for	General laboratory applications
DM Electrode/ Sensor	12 mm
measured size(s)	pH value temperature
pH MAX	14.0 pH
pH MIN	0.0 pH
Data transmission	Cable
Electrolyte	KCl
Shape lace	around
Material diaphragm	Ceramics
Material electrode / sensor	Borosilicate glass
Measuring unit	pH °C
Reference	Ag/AgCl
Type electrode / sensor	with cable connection
Type plug	BNC

