



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® pH electrode Flow T BNC NTC 30K

Special Price

€183.00 was

€235.00

**plus VAT &
Shipping**

Product Images



Description

The renowned Swiss manufacturer has developed a complete line of pH electrodes for glass-bodied laboratories. The new XS sensors are available with S7 head or cable fixation with BNC or DIN connector. A special polymer electrolyte/KCl or KCl electrolyte solution combined with a high performance reference system produces a full range of electrodes. Blue inner buffer provides a clear view of the contact with the inner pH glass membrane. All XS sensors can be connected to commercially available pH meters.

This model have been developed to facilitate the daily work in GLP condition. The internal serial number and certificate are quality tested.

Suitable for: Viscous samples and low-ion media

With temperature sensor NTC 30 $K\Omega$

Additional Information

No.	4-0941
Manufacturer (Brand)	neoLab
EAN	4058072189570
Transport temperature	Room temperature
Temp MAX	80 °C
Temp MIN	-10 °C
Length	120 mm
suitable for	General laboratory applications
DM Electrode/ Sensor	12 mm
Data transmission	Cable
Electrolyte	KCl
Shape lace	around
Material diaphragm	Polytetrafluoroethylene (PTFE. Teflon®)
Material electrode / sensor	Borosilicate glass
Measuring unit	pH °C
Reference	Ag/AgCl
Type electrode / sensor	with cable connection
Type plug	BNC
measured size(s)	pH value
pH MAX	14.0 pH
pH MIN	0.0 pH

