



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® pH electrode Micro P S8

€205.00
plus VAT &
Shipping

Product Images



Description

The renowned Swiss manufacturer has developed a complete line of pH electrodes for glass-bodied laboratories. The new XS sensors are available with S7 head or cable fixation with BNC or DIN connector.

A special polymer electrolyte/KCl or KCl electrolyte solution combined with a high performance reference system produces a full range of electrodes.

Blue inner buffer provides a clear view of the contact with the inner pH glass membrane.

All XS sensors can be connected to commercially available pH meters.

Some models have been developed to facilitate daily work in GLP condition.

The internal serial number and certificate are quality tested.

Suitable for: Small sample volumes <2 ml with protein

Additional Information

| | |
|-----------------------------|---------------------------------|
| No. | 4-0946 |
| Manufacturer (Brand) | neoLab |
| EAN | 4058072189624 |
| Transport temperature | Room temperature |
| Temp MAX | 100 °C |
| Temp MIN | 0 °C |
| Length | 120 mm |
| suitable for | General laboratory applications |
| DM Electrode/ Sensor | 12 mm |
| measured size(s) | pH value |
| pH MAX | 14.0 pH |
| pH MIN | 0.0 pH |
| DM Lace | 3.0 mm |
| Data transmission | Cable |
| Electrolyte | KCl |
| Shape lace | around |
| Material diaphragm | Ceramics |
| Material electrode / sensor | Borosilicate glass |
| Measuring unit | pH |
| Reference | Ag/AgCl |
| Type electrode / sensor | with cable connection |
| Type plug | S7 |

