



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® pH electrode Micro SpecialS8

€224.00
plus VAT &
Shipping

Product Images



Description

The renowned Swiss manufacturer has developed a complete line of pH electrodes for glass-bodied laboratories. The new XS sensors are available with S7 head or cable fixation with BNC or DIN connector.

A special polymer electrolyte/KCl or KCl electrolyte solution combined with a high performance reference system produces a full range of electrodes.

Blue inner buffer provides a clear view of the contact with the inner pH glass membrane.

All XS sensors can be connected to commercially available pH meters.

Suitable for: Measurements in NMR tubes

Additional Information

No.	4-0947
Manufacturer (Brand)	neoLab
EAN	4058072189631
Transport temperature	Room temperature
Temp MAX	100 °C
Temp MIN	0 °C
Length	220 mm
suitable for	General laboratory applications
DM Electrode/ Sensor	12 mm
DM Lace	3.0 mm
Data transmission	Cable
Electrolyte	KCl
Shape lace	around
Material diaphragm	Ceramics
Material electrode / sensor	Borosilicate glass
Measuring unit	pH
Reference	Ag/AgCl
Type electrode / sensor	with cable connection
Type plug	S7
measured size(s)	pH value
pH MAX	14.0 pH
pH MIN	0.0 pH

