



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657




## neolab® pH electrode Zava trode S8

**€250.00**  
**plus VAT &**  
**Shipping**

## Product Images

---



## Description

---

The renowned Swiss manufacturer has developed a complete line of pH electrodes for glass-bodied laboratories.

The new XS sensors are available with S7 head or cable fixation with BNC or DIN connector.

A special polymer electrolyte combined with a high performance reference system produces a full range of electrodes.

Blue inner buffer provides a clear view of the contact with the inner pH glass membrane.

All XS sensors can be connected to commercially available pH meters.

Insertion electrode/Suitable for cheese

## Additional Information

---

No.	4-0956
Manufacturer (Brand)	neoLab
EAN	4058072189723
Transport temperature	Room temperature
Temp MAX	60 °C
Temp MIN	0 °C
Length	190 mm
suitable for	Food sector
DM Electrode/ Sensor	12 mm
measured size(s)	pH value
pH MAX	14.0 pH
pH MIN	0.0 pH
Data transmission	Cable
Electrolyte	KCl
Shape lace	pointed
Material diaphragm	Ceramics
Material electrode / sensor	Borosilicate glass
Measuring unit	pH
Reference	Ag/AgCl
Type electrode / sensor	with cable connection
Type plug	S7

