



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® pH electrode 2-Pore K with stainless steel blade S8

€347.00
plus VAT &
Shipping

Product Images



Description

The renowned Swiss manufacturer has developed a complete line of pH electrodes for glass-bodied laboratories.

The new XS sensors are available with S7 head or cable fixation with BNC or DIN connector.

A special polymer electrolyte combined with a high performance reference system produces a full range of electrodes.

Blue inner buffer provides a clear view of the contact with the inner pH glass membrane.

All XS sensors can be connected to commercially available pH meters.

With stainless steel blade above the electrode tip

Additional Information

No.	4-0959
Manufacturer (Brand)	neoLab
EAN	4058072189754
Transport temperature	Room temperature
Temp MAX	60 °C
Temp MIN	0 °C
Length	110 mm
suitable for	Food sector
DM Electrode/ Sensor	35 mm
measured size(s)	pH value
pH MAX	14.0 pH
pH MIN	0.0 pH
Data transmission	Cable
Electrolyte	Polymer
Shape lace	pointed
Material diaphragm	Ceramics
Material electrode / sensor	Borosilicate glass stainless steel
Measuring unit	pH
Reference	Ag/AgCl
Type electrode / sensor	with cable connection
Type plug	S7

