



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Julabo Pt100 external sensor, 300 x 6 mm,
stainless steel, PTFE-coated, -75...+250 °C, 3 m**

**€585.00
plus VAT &
Shipping**

Product Images



Description

Pt100 external sensor, 300 x 6 mm \varnothing , stainless steel / PTFE coated, 3 m connection cable with kink protection. Applicable in temperature range from -75 °C to 250 °C. Without handle.

Additional Information

No.	JN-0356
Manufacturer (Brand)	Julabo
EAN	4058072846589
Transport temperature	Room temperature
Net weight	0.30000 kg
Country of origin	Germany
Customs tariff number	90251900

