



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Julabo Pt100 precision reference sensor, 180 x 4 mm D., stainless steel, 1 m connection cable

**€2,308.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Precision reference probe with integrated measurement electronics and RS232 connection for direct use with all JULABO calibration circulators. Measuring range -40.00...+300 °C, dimensions 180 x 4 mm ø, 1 m connection cable, immersion depth 90-140 mm.

## Additional Information

---

No.	JN-0350
Manufacturer (Brand)	Julabo
EAN	4058072846527
Transport temperature	Room temperature
Net weight	0.15000 kg
Country of origin	Germany
Customs tariff number	90251900

