



neoLab® Air/Gas Sensor NiCr-Ni

Product Images



Description

For measurements in seconds on surfaces, in liquids, air and gases as well as soft plastic media.

- 2 measuring ranges with different resolution
- Operation with any 2-pole NiCrNi probe
- 13 mm large, easy to read LCD display
- low power consumption, battery life up to 2000 hrs.
- ABS plastic housing

Technical data:

Total measuring range:

-65°C to +1150°C

Measuring range 1:

-65°C to +199,9°C

Resolution 0,1°C

Measuring range 2:

-65°C to + 1150°C

Resolution 1°C

Accuracy:

-65°C to +199,9°C: +/- 0,05%

-65°C to +1150°C: +/- 0.1%

Probe connection: 2-pin standard flat plug, for all common NiCrNi probes

Dimensions: 106 x 67 x 30 mm (HxWxD)

Weight: approx. 135 g

Additional Information

No.	4-1261
Manufacturer (Brand)	neoLab
EAN	4058072021467
Temp MAX	600 °C
Temp MIN	-65 °C
suitable for	Thermometer
Figure	2
Type probe	Luft-/ Gasfühler ohne Kabel
measured size(s)	Temperatur
Measuring unit	°C
for model measuring device	Digital-Sekunden-Thermometer (4-1252 und 4-1260)

