



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Surface probe NiCr-Ni

**€108.70**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

For measurements in seconds on surfaces, in liquids, air and gases as well as soft plastic media.

- 2 measuring ranges with different resolution
- Operation with any 2-pole NiCrNi probe
- 13 mm large, easy to read LCD display
- low power consumption, battery life up to 2000 hrs.
- ABS plastic housing

### Technical data:

Total measuring range:

-65°C to +1150°C

Measuring range 1:

-65°C to +199,9°C

Resolution 0,1°C

Measuring range 2:

-65°C to + 1150°C

Resolution 1°C

Accuracy:

-65°C to +199,9°C: +/- 0,05%

-65°C to +1150°C: +/- 0.1%

Probe connection: 2-pin standard flat plug, for all common NiCrNi probes

Dimensions: 106 x 67 x 30 mm (HxWxD)

Weight: approx. 135 g

## Additional Information

---

No.	4-1262
Manufacturer (Brand)	neoLab
EAN	4058072021474
Transport temperature	Room temperature
Temp MAX	900 °C
Temp MIN	-65 °C
suitable for	Thermometer
Figure	3
Measuring unit	°C
Type probe	Surface sensor without cable
for model measuring device	Digital second thermometer (4-1252 and 4-1260)
measured size(s)	Temperature

