



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Testo Waterproof Immersion/Piercing Probe NTC

€82.00
plus VAT &
Shipping

Product Images



Description

NTC temperature sensor Waterproof Measuring accuracy of up to ± 0.2 °C Measuring range: -50 to +150 °C The waterproof immersion/plunge probe with NTC sensor detects the temperature in liquid and semi-solid media quickly and precisely.

Additional Information

No.	TO-0216
Manufacturer (Brand)	Testo
old neoLab article no.	680931212
EAN	4058072373566
Transport temperature	Room temperature
Type probe	Insertion sensor Submersible sensor
measured size(s)	Temperature
Measuring unit	°C °F

