



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Testo immersion measuring tip TE type K (-200 to +1000 °C), 2 m thermo cable for WLAN data logger 2-5713

€148.00
plus VAT &
Shipping

Product Images



Description

Temperature probe for the radio data logger system testo Saveris 2 Suitable for model 2-T3 Flexible, low-mass immersion measuring tip, TE type K, oval cable Measuring range: -200°C to +1000°C

Additional Information

No.	TO-0747
Manufacturer (Brand)	Testo
old neoLab article no.	2-5731
EAN	4066292355439
Transport temperature	Room temperature
Temp MAX	1000 °C
Temp MIN	-200 °C
Length	500 mm
suitable for	Liquids
Type probe	Immersion sensor
DM Electrode/ Sensor	0.2 mm
measured size(s)	Temperature
Measuring unit	°C
for model measuring device	Saveris 2-T3 (1-5712)

