



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Testo Panel Sensor TE Type K (-50 to +205 °C) for WLAN Data Logger 2-5713

€54.00
plus VAT &
Shipping



Description

Temperature probe for the radio data logger system testo Saveris 2 Suitable for model 2-T3 panel-mounting probe with stainless steel sleeve, TE type K Measuring range: -50°C to +205°C

Additional Information

No.	TO-0794
Manufacturer (Brand)	Testo
old neoLab article no.	2-5730
EAN	4066292355866
Transport temperature	Room temperature
Temp MAX	205 °C
Temp MIN	-40 °C
Length	40 mm
suitable for	Liquids
DM Electrode/ Sensor	6 mm
Measuring unit	°C
Type probe	Built-in sensor
for model measuring device	Saveris 2-T3 (1-5712)
measured size(s)	Temperature

